

Reusable Digital Forehead Thermometer Products

Instructions for Use

How to Read

1. Apply to center of forehead just above the eyebrows. Hold in place until the color stops changing, generally 15 seconds. **Read the indicator while it is still on the forehead.**
2. Green indicates the correct temperature. If green does not appear, temperature is between blue and tan.

Intended Use

- Measure human forehead surface temperature and monitor temperature trends

Precautions

- Environmental temperature might affect the performance of this product
- Clean after use

Liquid Crystal thermometers are to be used on the forehead area only and estimate the core body temperature by measuring the temperature of the forehead. It is an ideal way to check for changes in temperature. A Liquid Crystal Forehead Thermometer is not as accurate for determining the actual core body temperature as a clinical (oral / rectal) thermometer. Instead it is designed to be an unbreakable, quick reading device for temperature screening between uses of a clinical thermometer.

Forehead skin temperatures are lower than oral temperatures.

Accurate to $\pm 1.0^{\circ}\text{F}$ / $\pm 0.6^{\circ}\text{C}$, the readings on a Liquid Crystal Forehead Thermometer have been adjusted to approximate the difference between forehead and oral temperatures during normal to feverish conditions. Liquid Crystal Forehead Thermometer products indicate an approximate oral temperature. If you think that your oral temperature is different from that shown by a Liquid Crystal Forehead Thermometer, you should check by using a clinical thermometer. If fever is present, consult your physician for appropriate treatment.

For External use only

This product contains no latex

Storage / Handling:

Products should be stored out of direct sunlight, dampness, or heat.

Cleaning:

Wipe clean using soap and water or rubbing alcohol. Do not soak.

Conforms to ASTM E-1061-01



LCR Hallcrest, 1911 Pickwick Lane, Glenview IL 60026 USA



LCR Hallcrest Ltd. Dock Road, Connah's Quay, UK CH5 4DS

Notified Body: Det Norske Veritas **CE** 0434