



Briteline

Product Data Sheet

Product Description: Briteline Urine Temperature Indicators

Standard Range: 89 to 101°F or 31.5°C to 38.5°C

Construction:

Adhesive Tape	Adhesive: Single-sided adhesive paper label (Fasson S-246)
Covering Film	Polyester
Color Change Material	Non toxic mixture of liquid crystal esters
Accuracy	±1°F / ± 0.6°C
Shelf Life	12 months from invoice date when stored at room temperature (i.e. 21°C (70°F) & 50% relative humidity).

Application Data: Apply temperature indicating label to outside of urine specimen cup. Reversible analog line moves along the scale to show detailed temperature reading in °F or °C.

Function: Display of temperature, as indicated by a bright line, of urine specimen cup. Label is adjusted with a 1°F / .6 °C offset.

Use: Ideal for use in clinics, hospitals, and doctors' offices.