



FeverScan Reusable

Product Data Sheet

Product Description: Reusable Forehead Temperature Monitor

Range: 94.6-104°F

Construction:

Covering Film	Polyester
Color Change Material	Non toxic mixture of liquid crystal esters
Accuracy	± 1.0°F / 0.6°C
Shelf Life	12 months from invoice date when stored at room temperature & humidity (i.e. 21°C (70°F) & 50% relative humidity).

Application Data: Peel temperature-indicating label from backing paper. Apply to forehead. Read temperature after colors stop changing. Green is correct reading. If green is not visible, subtract 1° from highest lit event.

Function: Continuous display of temperature. Core adjusted with 4.0°F / 2.2°C offset.

Latex Free