



# FeverScan Disposable

## Product Data Sheet

---

**Product Description:** Disposable Forehead Temperature Monitor

---

**Range:** 84-106°F

---

### Construction:

Adhesive	Hypoallergenic, pressure sensitive acrylate Carrier: Transparent 3 mil polyethylene film Release Liner: Bleached Kraft-Glassine paper, silicone coated on both sides, 3.5 mils Operating temp: -1 to 82 °C (continuous use)
Covering Film	Polyester
Color Change Material	Non toxic mixture of liquid crystal esters
Accuracy	± 1.0°F / 0.6 °C
Shelf Life	12 months from invoice date when stored at room temperature & humidity (i.e. 21°C (70°F) & 50% relative humidity).

---

**Application Data:** Peel temperature-indicating label from backing paper. Apply to forehead. Read temperature after colors stop changing. Green is correct reading. If green is not visible, subtract 1° from highest lit event.

**Function:** Continuous display of temperature. Core adjusted with 4.0°F / 2.2°C offset.

**Latex Free**



# FeverScan Reusable

## Product Data Sheet

---

**Product Description:** Reusable Forehead Temperature Monitor

---

**Range:** 94.6-104°F

---

### Construction:

Covering Film	Polyester
Color Change Material	Non toxic mixture of liquid crystal esters
Accuracy	± 1.0°F / 0.6°C
Shelf Life	12 months from invoice date when stored at room temperature & humidity (i.e. 21°C (70°F) & 50% relative humidity).

---

**Application Data:** Peel temperature-indicating label from backing paper. Apply to forehead. Read temperature after colors stop changing. Green is correct reading. If green is not visible, subtract 1° from highest lit event.

**Function:** Continuous display of temperature. Core adjusted with 4.0°F / 2.2°C offset.

**Latex Free**